

CERVICAL FUSION




ASPIS

Anterior Cervical Plate

5 levels available
Sizes from 22 mm to 110 mm

Prelordosed plates

Low profile : 1,9 mm
Width : 18 mm



+/- 14° angulation
ceph/caudal



Temporary fixation holes

One locking device for two
variable screws




Large visualization
windows

INSTRUMENTS



IMPLANTS



Length 22, 24, 26, 28, 30, 32, 34 mm

Length 36, 38, 40, 42, 44, 46, 48, 50 mm

Length 52, 55, 58, 61, 64, 67 mm

Length 70, 74, 78, 82, 86 mm

Length 90, 95, 100, 105, 110 mm

Screws \varnothing 4 and 4,5
Length 12, 14, 16, 18, 20, 22



4 Office, Seafieldpark,
Booterstown, Blackrock,
Co Dublin, Ireland

email: info@amedtherapeutics.com
web: www.amedtherapeutics.com