


Charis

PEEK PLIF Cages

Knowledge as leader
Human as vision



■ Lumbar Plif Cages



Characteristics

- **Bulleted Anterior Profile**
facilitates cage insertion into the intervertebral space
- **Retaining Teeth and Biconvex Profile**
prevent implant migration
- **Large Grafting Area**
to optimize bone fusion

CE
0499


amed
THERAPEUTICS

PEEK PLIF Cages

Charis



PEEK PLIF Cages

Knowledge as leader
Human as vision



| Implants Item No | Description | Implants Item No | Description |
|------------------|--------------------------------------|------------------|---------------------------------------|
| • CHA2274 | Length 22mm x Height 7mm (4 Degrees) | • CHA22104 | Length 22mm x Height 10mm (4 Degrees) |
| • CHA2278 | Length 22mm x Height 7mm (8 Degrees) | • CHA22108 | Length 22mm x Height 10mm (8 Degrees) |
| • CHA2574 | Length 25mm x Height 7mm (4 Degrees) | • CHA25104 | Length 25mm x Height 10mm (4 Degrees) |
| • CHA2578 | Length 25mm x Height 7mm (8 Degrees) | • CHA25108 | Length 25mm x Height 10mm (8 Degrees) |
| • CHA2284 | Length 22mm x Height 8mm (4 Degrees) | • CHA22114 | Length 22mm x Height 11mm (4 Degrees) |
| • CHA2288 | Length 22mm x Height 8mm (8 Degrees) | • CHA22118 | Length 22mm x Height 11mm (8 Degrees) |
| • CHA2584 | Length 25mm x Height 8mm (4 Degrees) | • CHA25114 | Length 25mm x Height 11mm (4 Degrees) |
| • CHA2588 | Length 25mm x Height 8mm (8 Degrees) | • CHA25118 | Length 25mm x Height 11mm (8 Degrees) |
| • CHA2294 | Length 22mm x Height 9mm (4 Degrees) | • CHA22124 | Length 22mm x Height 12mm (4 Degrees) |
| • CHA2298 | Length 22mm x Height 9mm (8 Degrees) | • CHA22128 | Length 22mm x Height 12mm (8 Degrees) |
| • CHA2594 | Length 25mm x Height 9mm (4 Degrees) | • CHA25124 | Length 25mm x Height 12mm (4 Degrees) |
| • CHA2598 | Length 25mm x Height 9mm (8 Degrees) | • CHA25128 | Length 25mm x Height 12mm (8 Degrees) |

Instruments



Produced exclusively for:

